




## General Purpose and Specific Electrode Cleaning Solutions

Electrodes can become dirty from use and will produce inaccurate results even as they read correctly in a pH buffer. Hanna's cleaning solutions eliminate impurities and residues that are left on the electrode surface after taking measurements.

Hanna general purpose and specific electrode cleaning solutions ensure maximum efficiency and accuracy of your sensors when used for its designated application. Electrode cleaning is a fast and effective routine that should be performed on a regular basis as a preventative measure against using a dirty electrode and to ensure that the junction is not clogged.



Electrode 1 has been properly cleaned before calibration.  
Electrode 2 has not been properly cleaned.

LINE ID:@neonics



ติดต่อบริษัท นีโอนิกส์ จำกัด

Tel: 098-479-5684 หรือ 061-8268939

E-mail: sale@neonics.co.th หรือ sale@tools.in.th

## General Electrode Cleaning

Clean the sensing portion of your electrodes weekly to prevent fouling and to maintain accuracy. Immerse the electrode in the proper cleaning solution for at least 15 to 20 minutes and rehydrate in storage solution before use.

### General Use Electrode Cleaning Solutions

| Code       | Application               | Package           |
|------------|---------------------------|-------------------|
| HI70000P   | rinsing                   | 20 mL sachet (25) |
| HI7061L    | general purpose           | 500 mL bottle     |
| HI7061-023 | general purpose (GroLine) | 230 mL bottle     |
| HI7061-012 | general purpose (GroLine) | 120 mL bottle     |
| HI7073L    | proteins                  | 500 mL bottle     |
| HI7073M    | proteins                  | 230 mL bottle     |
| HI7074L    | inorganic substances      | 500 mL bottle     |
| HI7074M    | inorganic substances      | 230 mL bottle     |
| HI7077L    | oil and fats              | 500 mL bottle     |
| HI7077M    | oil and fats              | 230 mL bottle     |
| HI8061L    | general purpose           | 500 mL FDA bottle |
| HI8073L    | proteins                  | 500 mL FDA bottle |
| HI8077L    | oil and fats              | 500 mL FDA bottle |



### Specific Cleaning for a Top Performing Sensor

In many applications, electrodes become contaminated from use and produce inaccurate results. Since these contaminants cannot be removed during normal use, special cleaning solutions are needed.

The Cleaning Series ensures maximum efficiency and accuracy of your sensors when used for its designated application

### Specific Cleaning Solutions - Bottles

| Code     | Description   | Size   |
|----------|---|--------|
| HI70630L | acid cleaning solution for meat grease and fats (food industry)                         | 500 mL |
| HI70631L | alkaline cleaning solution for meat grease and fats (food industry)                     | 500 mL |
| HI70632L | cleaning and disinfection solution for blood products                                   | 500 mL |
| HI70635L | cleaning solution for wine deposits (winemaking)  | 500 mL |
| HI70636L | cleaning solution for wine stains (winemaking)  | 500 mL |
| HI70640L | cleaning solution for milk deposits (food industry)                                     | 500 mL |
| HI70641L | cleaning and disinfection solution for dairy products (food industry)                   | 500 mL |
| HI70642L | cleaning solution for cheese residues (food industry)                                   | 500 mL |
| HI70643L | cleaning and disinfection solution for yogurt products (food industry)                  | 500 mL |
| HI70663L | cleaning solution for soil deposits (agriculture)                                       | 500 mL |
| HI70664L | cleaning solution for humus deposits (agriculture)                                      | 500 mL |
| HI70670L | cleaning solution for salt deposits (industrial processes)                              | 500 mL |
| HI70671L | cleaning and disinfection solution for algae, fungi and bacteria (industrial processes) | 500 mL |
| HI70681L | cleaning solution for ink stains  | 500 mL |

### Specific Cleaning Solutions - Sachets

| Code      | Description  | Qty/Size   |
|-----------|--|------------|
| HI700601P | general purpose cleaning solution for laboratories                     | 20 mL (25) |
| HI700630P | acid cleaning solution for meat grease and fats (food industry)        | 20 mL (25) |
| HI700635P | cleaning solution for wine deposits (winemaking)                       | 20 mL (25) |
| HI700636P | cleaning solution for wine stains (winemaking)                         | 20 mL (25) |
| HI700640P | cleaning solution for milk deposits (food industry)                    | 20 mL (25) |
| HI700641P | cleaning and disinfection solution for dairy products (food industry)  | 20 mL (25) |
| HI700642P | cleaning solution for cheese residues (food industry)                  | 20 mL (25) |
| HI700643P | cleaning and disinfection solution for yogurt products (food industry) | 20 mL (25) |
| HI700661P | general purpose cleaning solution for agriculture                      | 20 mL (25) |
| HI700663P | cleaning solution for soil deposits (agriculture)                      | 20 mL (25) |
| HI700664P | cleaning solution for humus deposits (agriculture)                     | 20 mL (25) |
| HI700670P | cleaning solution for salt deposits (industrial processes)             | 20 mL (25) |

LINE ID:@neonics



ติดต่อบริษัท นีโอนิกส์ จำกัด

Tel: 098-479-5684 หรือ 061-8268939

E-mail: sale@neonics.co.th หรือ sale@tools.in.th